



User Guide

Immobilon[®]-P Transfer Membrane

For High Sensitivity Immunodetection

Notice

The information in this document is subject to change without notice and should not be construed as a commitment by Merck Millipore Ltd. ("Millipore") or an affiliate. Neither Merck Millipore Ltd. nor any of its affiliates assumes responsibility for any errors that may appear in this document.

The M logo, Millipore, Immobilon, bløk, SNAP i.d., and Milli-Q are registered trademarks of Merck KGaA, Darmstadt, Germany. Luminata is a trademark of Merck KGaA.

All trademarks of third parties are the property of their respective owners.

© 2015 EMD Millipore Corporation. Billerica, MA, U.S.A. All rights reserved.

PR04715, Rev. 04/15

Contents

Introduction.....	1
Guidelines for Working with Immobilon®-P membrane.....	2
Materials Recommended for Western Blotting	2
Protein Transfer.....	3
Immunodetection.....	4
Common protocols used in Western blotting and Immunodetection.....	4
Membrane Wetting	4
Semi-dry Transfer.....	4
Protein Visualization (Optional).....	5
Immunodetection.....	5
SNAP i.d.® 2.0 Immunodetection Using Vacuum Filtration	6
Guidelines for Choosing an Immobilon® PVDF Membrane	6
Ordering Information	7
Technical Assistance.....	8
Standard Warranty	8

Introduction

Immobilon®-P transfer membrane is a polyvinylidene fluoride (PVDF) microporous membrane used for transfer of proteins from a variety of gel matrices. When compared to a nitrocellulose membrane, it has improved handling characteristics and staining capabilities, increased solvent resistance, and a higher signal-to-noise ratio for enhanced sensitivities.

This hydrophobic membrane has a nominal pore size of 0.45 micron (µm) and is capable of binding a wide range of molecular weight proteins. While Immobilon®-P membrane is capable of binding proteins smaller than 20,000 in molecular weight, Immobilon®-P⁵⁰ 0.2 micron (µm) membrane could be considered for proteins smaller than 20,000, because it has more surface area and higher binding capacity for small proteins.

Immobilon®-P membrane has excellent protein retention, high physical strength, and broad chemical compatibility, making it ideal for a variety of staining applications and reprobing in immunodetection methods. This user guide describes some of the most common techniques for Western blotting. For more detailed information please refer to [TPO01EN](#), the "Protein Blotting Handbook" (available at www.millipore.com).

Table 1. Immobilon-P® Membrane Properties and Applications

Composition	PVDF	
Pore size	0.45 µm	
Phobicity	Hydrophobic	
Protein binding capacity	Insulin: 160 µg/cm ² Bovine serum albumin (BSA): 215 µg/cm ² Goat IgG: 294 µg/cm ²	
Applications	Binding assays Dot/slot blotting Glycoprotein visualization Lipopolysaccharide analysis Mass spectrometry Amino acid analysis N- terminal protein sequencing	
Detection methods*	Chemiluminescent (Luminata™ and Immobilon® HRP substrates) Chromogenic (Calbiochem® TMB) Radioactive	
Protein visualization methods	Transillumination	
	Reversible stains	Ponceau-S Copper phthalocyanine tetrasulfonic acid, tetrasodium salt (CPTS) Toluidine blue Sypro® blot stains
	Irreversible stains	Coomassie™ Brilliant Blue dye Amido black India ink Colloidal gold
* For fluorescence detection methods, low-autofluorescent Immobilon®-FL membrane is recommended.		

Guidelines for Working with Immobilon®-P membrane

- Always wear gloves when handling the membrane, in order to avoid fingerprints.
- Use blunt forceps to prevent membrane damage.
- Keep the patapar (blue paper) with the membrane during cutting or handling, but discard when wetting the membrane.
- Handle with care to avoid scratches on the membrane surface. Do not fold the membrane.
- Hydrophobic Immobilon®-P membrane must be wet in an alcohol solution (> 50% v/v methanol, ethanol, or isopropanol) before use. Once the membrane is wet, it changes from opaque to semi-transparent.
- After protein transfer, wash the blot with Milli-Q® water to eliminate any gel residues.
- Blots can be air dried and stored at 4° C for several months (for later use) or they can be used immediately.


Materials Recommended for Western Blotting

- Immobilon®-P membrane cut to the dimensions of the gel
- Alcohol (> 50% methanol, ethanol, or isopropanol) for wetting dry membrane
- Milli-Q® water
- Transfer buffer: 25 mM Tris-base, 192 mM glycine, pH 8.3, 10% alcohol for tank transfer or 48 mM Tris, 39 mM glycine, pH 9.2, 10% alcohol for semi-dry transfer)
- Sheets of filter paper, cut to the dimensions of the gel and soaked in transfer buffer for at least 30 seconds
- Blocking buffer: bløk®-CH buffer (cat. no. WBAVDCH01) or 0.5–5% (w/v) blocking agent (bovine serum albumin, casein, nonfat dry milk) in wash buffer
- Wash buffer: Phosphate-buffered saline (PBS) or Tris-buffered saline (TBS) containing 0.5–0.1% Tween®-20 surfactant (PBST or TBST)
 - PBS: 10 mM sodium phosphate, pH 7.2, 0.9% NaCl
 - TBS: 10 mM Tris, pH 7.4, 0.9% NaCl
- Primary antibody (specific for the protein of interest), diluted in blocking buffer or wash buffer
- Secondary antibody (specific for the primary antibody), labeled with a detection enzyme (e.g., horseradish peroxidase [HRP] or alkaline phosphatase [AP]), diluted in blocking buffer or wash buffer

Protein Transfer

Proteins can be transferred to Immobilon®-P membrane by two common electro-transfer methods: tank and semi-dry transfer. Table 2 describes the general conditions and major differences for the two methods.

Table 2. Transfer Methods



Conditions	Membrane gel/ stack immersed in buffer	Filter paper soaked in buffer
Buffer volume	0.5 L or higher depending on system	~ 0.05 L per mini-gel
Transfer time	Slow (1 or more hours)	Fast (7–45 minutes)
Typical run	Constant voltage	Constant current
Gel equilibration time	Not necessary but recommended	Minimum of 15 minutes equilibration per mini-gel
Required buffer system	Continuous (single buffer)	Continuous (single buffer) or Discontinuous (3 buffers)
Typical continuous buffer name/composition	Towbin buffer, pH 8.3 25 mM Tris, 192 mM glycine	Bjerrum–Shafer–Nielsen buffer, pH 9.2 48 mM Tris, 39 mM glycine
Typical discontinuous buffer name/composition	N/A	Anode buffer I: 300 mM Tris, pH 10.4 Anode buffer II: 25 mM Tris, pH 10.4 Cathode buffer: 25 mM Tris, 40 mM amino-caproic acid, pH 9.4
% Alcohol*	10–20% alcohol	10–20% alcohol
% SDS*	0.02–0.05% SDS	0.025–0.05% SDS

* Alcohol (methanol, ethanol, or isopropanol) in the transfer buffer has two important functions; it stabilizes the gel dimensions and strips complexed sodium dodecyl sulfate (SDS) from protein molecules, improving protein binding to the membrane. However, for large proteins, or proteins that exhibit solubility problems, it is recommended that the alcohol concentration be decreased and that a small amount of SDS be added to the transfer buffer. This improves protein elution from the gel while maintaining protein solubility during the transfer process.

Immunodetection

Immunodetection is an antibody-based method that allows the detection, identification, and quantitation of a protein or antigen in the blotting membrane. The general protocol follows these six general steps:

1. Block unoccupied membrane sites to prevent nonspecific binding of antibodies.
2. Incubate the membrane with a primary antibody that binds to the protein of interest.
3. Wash to remove any unbound primary antibody.
4. Incubate the membrane with a conjugated secondary antibody, which binds to the first antibody.
5. Wash to remove any unbound secondary antibody.
6. Incubate the membrane with a substrate that reacts with the conjugated secondary antibody to reveal the location of the protein.

Standard immunodetection takes at least 4 hours and is widely used, but the SNAP i.d.[®] Protein Detection System (www.millipore.com/snapwb) can perform the same process with significant time savings (Table 3).

Table 3. Comparison of Standard vs. SNAP i.d.[®] Immunodetection

	Standard Immunodetection	SNAP i.d. [®] Immunodetection
1. Block membrane	1 hour	10 seconds
2. Incubate with primary antibody	1 hour	10 minutes
3. Wash membrane	3 times, 5 minutes each	4 times, ~10 seconds each
4. Incubate with secondary antibody	1 hour	10 minutes
5. Wash membrane	3 times, 5 minutes each	4 times, ~10 seconds each
Total time	4 hours	22 minutes

Common protocols used in Western blotting and Immunodetection

Membrane Wetting

1. Wet the dry membrane in alcohol (>50% methanol, ethanol, or isopropanol) for 10–20 seconds, or until it changes from a opaque white to uniform, translucent gray.
2. Immerse the membrane in Milli-Q[®] water for 1–2 minutes to displace the alcohol.
3. Equilibrate the membrane in transfer buffer for 2–3 minutes or until ready to use.

CAUTION: Once the membrane has been wet out, do not allow it to dry out. It can be kept in buffer until protein transfer. If the membrane dries out (turns opaque white) even partially, it must be wet out again (steps 1–3).

Semi-dry Transfer

1. Resolve the protein mixture on a 1D or 2D polyacrylamide gel.
2. Immerse the gel in the transfer buffer and allow it to equilibrate for 10–15 minutes.
3. Assemble the transfer stack according to manufacturer's instructions for the transfer apparatus used.

CAUTION: To ensure an even transfer, remove air bubbles by carefully rolling a clean pipette or blot roller over the surface of each layer in the stack. Do not apply excessive pressure, as this may damage the gel and membrane.

Semi-dry Transfer, continued

4. Transfer proteins according to transfer apparatus manufacturer's instructions.
5. Remove the blot from the transfer system and rinse the membrane briefly in Milli-Q® water to remove gel debris. The blot may be air dried for storage, or it can be use immediately for the immunodetection step.

NOTE: Drying the blot before immunodetection may enhance the binding of some proteins and reduce background noise.

Protein Visualization (Optional)

To visualize the protein transfer efficiency, Immobilon®-P membrane may be stained with any reversible blot stain compatible with immunodetection (e.g., Ponceau-S, toluidine blue, CPTS, Sypro® blot stains) or viewed by transillumination using a light box. For a list of reversible and nonreversible compatible stains and protocols refer to [TP001EN](#), the "Protein Blotting Handbook" at www.millipore.com.

Immunodetection

The following is a general protocol for immunodetection with Immobilon®-P membrane. Some of the critical factors for obtaining a "perfect" Western blot such as protein concentration, blocking solution, and antibody concentration may require optimization.

Standard Immunodetection

1. If blot has been dried, rewet it in alcohol (> 50% methanol, ethanol, or isopropanol) for 15 seconds or until it changes from opaque white to translucent gray.
2. Rinse the blot in Milli-Q® water for 1 minute.
3. Place the blot in blocking buffer and incubate for 1 hour with gentle agitation. Dilute the primary antibody in wash or blocking buffer.
4. Place the blot in diluted primary antibody solution and incubate for 1 hour with gentle agitation.
5. Wash the blot with wash buffer (tris- or phosphate-buffered saline solution, supplemented with Tween® 20 surfactant (TBST or PBST)) 3–5 times for 5 minutes each. Prepare secondary antibody in wash or blocking buffer.
6. Place the blot in diluted enzyme-labeled secondary antibody solution and incubate for 1 hour with gentle agitation.
7. Wash the blot with wash buffer 3–5 times for 5 minutes each.
8. Place the blot into a clean container and add the appropriate detection reagent (HRP, AP, or chromogenic).
9. Incubate 1–5 minutes, according to the detection reagent manufacturer's instructions.
10. For HRP or AP chemiluminescent reagents, expose blot to x-ray film or acquire the image using a digital imaging system. For chromogenic detection, add the reagent and wait until signal appears.

SNAP i.d.[®] 2.0 Immunodetection Using Vacuum Filtration

(refer to the [SNAP i.d.[®] Protein Detection System User Guide](#) for full protocol details)

1. If blot has been dried, rewet it in > 50% alcohol (methanol, ethanol, or isopropanol) for 15 seconds, or until it changes from opaque white to translucent gray. Prepare all the required solutions and antibodies ahead of time.

NOTE: Antibodies should be 3 to 5 times more concentrated than in standard immunodetection, but in volumes of 2.5 to 10 mL depending of the blot size/blot holder.

2. Wet the SNAP i.d.[®] blot holder in Milli-Q[®] water and assemble the blot with the protein side down.
3. Using the blot roller, remove all air bubbles and excess water, and insert the blot holder inside the SNAP i.d.[®] frame.
4. Block by adding 15–30 mL of blocking solution and immediately turn the vacuum on.
5. Depending of the blot holder size used, add 2.5 to 10 mL of diluted primary antibody and incubate for 10 minutes.
6. Turn vacuum on to flush the antibody, then with the vacuum still on, wash 4 times with 15–30 mL of wash buffer.
7. Turn vacuum off and add 2.5–10 mL (depending of the blot holder size used) of diluted secondary antibody and incubate for 10 minutes.
8. Turn vacuum on to flush the antibody, then wash 4 times with 15–30 mL of wash buffer.
9. Remove the blot from the blot holder and continue with the detection method of choice (chemiluminescence or chromogenic).

Guidelines for Choosing an Immobilon[®] PVDF Membrane

The following table provides general guidelines for choosing the appropriate membrane for a specific post-Western blot application. Due to variations in protein properties such as charge density, conformation, and hydrophobicity, not all proteins behave the same way on a given membrane surface. Experiments with a variety of Immobilon[®] membranes may be necessary to optimize results for your specific application.

Application after Western blotting	Membrane of choice for most proteins
General immunodetection	Immobilon [®] -P
Amino acid analysis	Immobilon [®] -P
Immunodetection of low molecular weight or low-abundance proteins	Immobilon [®] -P ^{SO}
Sequencing of low molecular weight or low-abundance proteins	Immobilon [®] -P ^{SO}
Fluorescence immunodetection and chemifluorescence methods	Immobilon [®] -FL

Ordering Information

See the Technical Assistance section for contact information. You can purchase these products on-line at www.millipore.com/products.

Immobilon®-P Membrane (0.45 µm pore size) for general Western blotting applications

Size	Qty/Pk	Catalogue Number
26.5 × 375 cm roll	1	IPVH00010
26 × 26 cm sheet	10	IPVH304F0
20 × 20 cm sheet	10	IPVH20200
15 × 15 cm sheet	10	IPVH15150
10 × 10 cm sheet	10	IPVH10100
9 × 12 cm sheet	10	IPVH09120
8.5 × 13.5 cm sheet	10	IPVH08130
8 × 10 cm sheet	10	IPVH08100
7 × 8.4 cm sheet	50	IPVH07850

Immobilon®-P^{SO} Membrane (0.2 µm pore size) for blotting applications of proteins with molecular weights less than 20,000

Size	Qty/Pk	Catalogue Number
26.5 × 375 cm roll	1	ISEQ00010
26 × 26 cm sheet	10	ISEQ26260
20 × 20 cm sheet	10	ISEQ20200
15 × 15 cm sheet	10	ISEQ15150
10 × 10 cm sheet	10	ISEQ10100
9 × 12 cm sheet	10	ISEQ09120
8.5 × 13.5 cm sheet	10	ISEQ08130
8 × 10 cm sheet	10	ISEQ08100
7 × 8.4 cm sheet	50	ISEQ07850

Immobilon®-FL Membrane (0.45 µm pore size) for fluorescence detection applications

Size	Qty/Pk	Catalogue Number
26 cm × 375 cm roll	1	IPFL00010
20 × 20 cm sheet	10	IPFL20200
10 × 10 cm sheet	10	IPFL10100
7 × 8.4 cm sheet	10	IPFL07810

Technical Assistance

For more information, contact the office nearest you. In the U.S., call 1-800-221-1975. Outside the U.S., go to our web site at www.millipore.com/offices for up-to-date worldwide contact information. You can also visit the tech service page on our web site at www.millipore.com/techservice.

Standard Warranty

The applicable warranty for the products listed in this publication may be found at www.millipore.com/terms ("Conditions of Sale").